




Caprosyn™ Sutures



Tensile strength:


- Caprosyn™ retains 50-60% knot strength at 5 days post implantation and a minimum of 20-30% of knot strength at 10 days post implantation. All tensile strength is essentially lost by 21 days.

Rapid absorption:

- Caprosyn™ synthetic absorbable suture is completely absorbed by 56 days.

Use:

- Soft tissue closure where prolonged support is not needed (eg. surgery of the oral cavity)



SC691
SC692
GC122
GC123



Polyglytone™ 6211 the fastest absorbing monofilament synthetic suture

Main characteristics:

- Strong and secure wound approximation during critical wound healing period
- Rapid absorption
- Excellent handling
- Soft and flexible texture
- Excellent knot security