



Skin Staplers

Royal™ skin stapler

The Royal™ skin stapler is a high quality stapler with a convenient hand grip to allow a clear line of sight for accurate positioning of staples. The stapler is presented in 2 sizes each with 35 staples (35 F and 35 W).

Main characteristics:

- Clear line of sight for accurate positioning of staples
- Quick wound closure
- Complete closure of staple to minimise turning
- Aesthetic and non-reactive

Royal™ 35 (standard)	Royal™ 35 W (wide)
35 Staples for thin skin	35 Staples for thick skin
	

Appose™ skin stapler

The Appose™ skin stapler also comes in two sizes, 35 R and 35 W, each with 35 staples.



Appose™ 35 (standard)	Appose™ 35 W (wide)
35 Staples for thin skin	35 Staples for thick skin
	

Multifire Premium™ skin stapler

The Multifire Premium™ skin stapler is a high quality single use re-loadable skin stapler. Each loading unit contains 35 staples and the instrument can be re-loaded twice, giving a total of 105 staples (3 cartridges). Staples come in two sizes: normal and wide.

Main characteristics:

- Quick wound closure
- Re-loadable skin stapler (giving a total of 105 staples)
- 360 degrees rotating head, maximising access to wound edges
- Stainless steel staples
- Very suitable for equine surgery

Multifire Premium™ 35 (standard)	Multifire Premium™ 35 W (wide)
35 Staples for thin skin	35 Staples for thick skin
	



AGRAF002
AGRAF001



Skin staplers provide a very efficient method of skin closure requiring minimal chemical restraint and analgesia to insert and remove. They are ideal for quick procedures such as cut pads and replacing sutures removed by an animal's teeth. Skin staples are also ideal if you need to shorten anaesthetic time as they are faster to insert compared to suture material.

Skin staple extractor

The skin staple extractor is a high quality staple extractor, facilitating pain-free staple removal. It can be used with Royal™, Appose™ and Multifire Premium™ staples.



OTE002



Skin Staple Removers