

Thoracoabdominal instruments are useful in a wide variety of thoracic and abdominal procedures. TA™ devices can be totally disposable, “multifire”, or come in the form of a surgical stainless steel model that can be steam sterilised.

**Main applications in veterinary surgery:**

- Partial and complete lung and liver lobectomy
- Resection of tumours
- Gastro-intestinal resection and anastomosis
- Enterotomy closure

**Reusable TA™ stainless steel staplers**

Stapler	Disposable Loading Unit Colour	Length of staple line (mm)	Number of staple rows	Staple size (mm)
TA™30V3	White	30	3	2.0
TA™30	Blue	30	2	3.5
	Green	30	2	4.8
TA™55	Blue	55	2	3.5
	Green	55	2	4.8
TA™90	Blue	90	2	3.5
	Green	90	2	4.8

These devices accept a one-piece disposable loading cartridge, that contains the staple unit, the opposing anvil and retaining pin. The function of the TA™ stapler is to place double or triple staggered rows of “B” shaped staples.



TY030470  
TY030715  
TY3975



**Single use reloadable DST Series™ TA™ staplers**

DST Series™ TA™ staplers have an ergonomic, easy to use design. Each DST Series™ TA™ stapler can be re-loaded seven times for a total of eight firings, as each stapler is initially packed with a loaded and ready-to-use cartridge. These staplers contain a dual action handle for closing and firing staples, with a pre-clamp position that allows the veterinary surgeon to properly position and check the tissue prior to firing the stapler.

Stapler	Order code	Length of Staple Line (mm)	Loading Unit	Order Code	Disposable Loading Unit Colour	Staple size (mm)	Number of staple rows
<b>Autosuture TA™ DST Series™ Single use re-loadable stapler</b>							
TA30V3S	408877	30	TA30V3L	408876	White	2.5	3
TA3035S	408873	30	TA3035L	0808511	Blue	3.5	2
TA3048S	408875	30	TA3048L	408874	Green	4.8	2
TA4535S	408879	45	TA4535L	408878	Blue	3.5	2
TA4548S	408881	45	TA4548L	408880	Green	4.8	2
TA6035S	0808819	60	TA6035L	0808478	Blue	3.5	2
TA6048S	0808830	60	TA6048L	408882	Green	4.8	2
TA9035S	408884	90	TA9035L	408883	Blue	3.5	2
TA9048S	408885	90	TA9048L	0808456	Green	4.8	2